

For CHO cells  
Recombinant Protein / Virus Production Tool

# Gxpress CHO Transfection & Medium Kit

## ■ Gxpress CHO Transfection & Medium Kit

Package Contents (Catalog No. GXCHO-MAK-0010 ¥72,900)

- Gxpress CHO TF Reagent 1.2 mL x 2
- CHO Complexing Solution 40 mL x 2
- Gxpress CHO Enhancer 120 mL x 1
- CH400AZ medium 1000 mL x 1




### Description


Mix gene and reagent, and add to culture medium. Next day, enhancer reagent can be added to produce desired recombinant protein in a short period of time.

## Gxpress CHO Transfection & Medium Kit


### 1) Growth assay of CHO cells in shaker culture.



### Protocol



### 2) Short-term production of human IgG using Gxpress CHO Transfection & Medium Kit.



Product Name	Catalog Number	Enough to Transfect	Price
Gxpress CHO Transfection & Medium Kit *1	GXCHO-MAK-0010	1 L culture	¥72,900
Gxpress CHO Transfection & Medium Kit II *2	GXCHO-MK-0010	1 L culture	¥71,100
Gxpress CHO Transfection Kit	GXCHO-RK-0010	1 L culture	¥53,500
Gxpress CHO TF Reagent	GXCHO-TF-0010	1 L culture	¥42,900
Gxpress CHO Enhancer	GXCHO-EN-0010	1 L culture	¥16,500

\*1 It is a ready-to-use product including CH400AZ medium (with Ala-Gln) to which L-alanyl-L-glutamine is added as a product component.

\*2 It is a product including CH400 medium (w/o Ala-Gln or Gln) that can control glutamine concentration during culture.

Gmep Incorporated

Hyakunen Park, Kurume, Fukuoka, Japan 8390864  
WEB : <https://gmep.co.jp> TEL : +81-0942-65-4463

○ The listed price is as of October 2020. Please note in advance that the prices and product configurations described may change.  
○ These items are for research use only. Please do not use for human and animal diagnosis, treatment, other conditions.